



DUO-JET[®] II

LONG-TERM HEMODIALYSIS CATHETER



DUO-JET[®] II

LONG-TERM HEMODIALYSIS CATHETER



10F CATHETER SET

DJLT2000L	10F x 52CM DUO-JET [®] II SINGLE LUMEN CATHETER, LEFT SIDE APPROACH	5/BOX
DJLT2000R	10F x 52CM DUO-JET [®] II SINGLE LUMEN CATHETER, RIGHT SIDE APPROACH	5/BOX

Catheter Set Contents: (2) Catheter (2) Introducer Needle (2) Guidewire with Advancer (2) Tunneler (2) Peelable Introducer (1) Coaxial Introducer Sheath (1) Arterial Extension Set (1) Venous Extension Set (2) Clamp (2) Catheter Plug (2) End Cap (1) Catheter Securement Device (1) Patient ID Card (1) Patient Information Packet

IMPORTANT RISK INFORMATION

Indications For Use: The Duo-Jet[®] II Catheter is indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis.

Contraindications: Known or suspected allergies to any of the components of the catheter or the kit. This device is contraindicated for patients exhibiting severe, uncontrolled coagulopathy or thrombocytopenia.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications. For patient information please visit www.medcomp.net/patientinformation



Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201 F: 215.256.1787
www.medcomp.net

ISO13485
Certified Quality System

PN2674BSI Rev. A 10/22